



Solo™

BENSON MEDICAL INSTRUMENTS

Occupational Database Management System

Solo™ gives you the tools to keep your hearing conservation program in compliance with OSHA regulations. You can produce reports and notification letters in seven languages based on subjects individual results.

Ease of use is evident in Solo's™ design. You'll find all your subjects and audiograms conveniently displayed on the main program window. Search for subjects by name or Social Security number. Click on an audiogram to see results compared to the baseline audiogram.

Transferring an audiogram from your audiometer has never been easier. Press a button to automatically transfer the data and print the report of your choice, even while you're editing data.

The selection of reports is comprehensive. Display test results several different ways with the many Subject reports and letters. Company reports show results for the whole company. See how Solo™ will help you administer an effective hearing conservation program.

- Creates notification letters in seven languages
- Produces individual and summary reports
- Manages programs for multiple companies
- Calculates impairments, standard threshold shifts, recordable shifts, and medical referrals
- Meets OSHA requirements—easily and economically
- Customizes data look-up tables
- Transfers hearing data from your audiometer

PO Box 190 ~ 14221 Helen Lane
Mountain WI 54149-0190
(800) 356-5239 Voice ~ (800) 644-7408 Fax
srmicom Website ~ info@srmicom